



# Simplify with SimplyGo

Now one POC meets the portable oxygen needs of nearly all of your patients

**PHILIPS**  
RESPIRONICS



# Treating your oxygen patients has never been easier

SimplyGo is the only portable oxygen concentrator (POC) to offer continuous flow and pulse-dose delivery in a single device weighing only 10 pounds.

That means you can now manage the therapy and lifestyle needs of nearly all of your patients with just one POC. Designed to be rugged and easy to use, SimplyGo is truly a new breed of lightweight and attractive portable concentrator. Simplify patient management – and your business – with the SimplyGo portable oxygen concentrator from Philips Respironics.

A new breed of portable concentrator that simplifies the care and management of nearly all oxygen patients.





# SimplyGo simply...



## Meets most patients' needs

SimplyGo packs more of what you need in one POC, without packing on unwanted pounds. In fact, SimplyGo is the only continuous flow POC to weigh in at ten pounds. And for patients who need higher pulse settings at night, the Sleep mode provides increased trigger sensitivity and a softer pulse for exceptional comfort. There's simply no POC more flexible than SimplyGo.

## Redefines ruggedness

When it comes to toughness, SimplyGo is no pushover. With a long-life compressor, high-quality parts, a high-impact resistant design, and oversized cart wheels, this POC was built to set new standards of ruggedness. During development, SimplyGo was rigorously tested and subjected to extreme conditions including impacts, vibrations, and temperatures to deliver reliable performance day in and day out.


## Fits your business

SimplyGo's warranty and service programs are all about flexibility. With SimplyGo, you have the choice of a two-, three-, or five-year warranty. Pick the plan that's best for you. While some POCs must be returned to the manufacturer for service, SimplyGo's easily accessible parts and comprehensive service tools allow your technicians to perform the work right in your shop if you prefer.

# SimplyGo for today and tomorrow

SimplyGo has more than twice the oxygen output of any POC weighing ten pounds or less. That means that SimplyGo is able to keep up with your patients' portable oxygen needs today and tomorrow.

## Oxygen output of POCs weighing 10 pounds or less



## Powered for success

One of the keys to SimplyGo's portability is its lightweight and compact battery. SimplyGo's battery is immediately accessible and can be easily changed by your patients. There are no awkward external battery belts or connectors that get in the way, or hard-to-maintain internal batteries.

### SimplyGo's battery duration

| Pulse dose setting of 2* | Continuous flow 2 LPM |
|--------------------------|-----------------------|
| 3.5 hours                | .7 hours              |

\*Based on 20 breaths per minute.

Actual duration results may vary.



The basic SimplyGo system comes with one Lithium Ion battery so you don't have to buy more batteries than you need. But if your patients want more power, extra batteries fit easily into the zippered pouch or handy accessory bag. Patients want more time? Simply take along more batteries.

## Intended to streamline your business...

Since SimplyGo meets the needs of nearly all of your patients, you can greatly simplify your inventory management and use less shelf space, helping you to lower costs and streamline patient services. And, SimplyGo's ability to provide continuous flow in a small, lightweight package is appealing to referral sources and patients.



## ...and be a team player

When SimplyGo is teamed up with EverFlo – our low maintenance, lightweight, stationary oxygen concentrator for use in the home – you can create an oxygen-generating system that will meet the needs of nearly all of your patients while helping to reduce your cost of providing oxygen therapy.

Our SimplyGo/EverFlo duo helps:

- Eliminate recurring oxygen delivery costs and cylinder hassles
- Reduce maintenance costs
- Ease inventory management
- Minimize patient calls



SimplyGo is built with pride at our award-winning facility in Kennesaw, Georgia.



## One attractive package

SimplyGo's non-clinical style and ease-of-use make it appealing to patients. An attractive and highly functional carrying case allows users to either carry SimplyGo over the shoulder or place it on the easy-to-pull, easy-to-store cart. Power from a Lithium Ion battery or standard AC/DC power cords add convenience and flexibility. So, whether it's a walk across the street, a ride across town, or a trip across the country, we've made it easier than ever for you to help patients – and yourself – simplify with SimplyGo.



## The SimplyGo portable oxygen concentrator

Comes with everything that you need including a carrying case, mobile cart, one rechargeable Lithium Ion battery, DC and AC power cords, accessory bag and user manual.

- Attractive, sleek design
- Easy-to-remove battery
- Intuitive, easy-to-read screen
- Convenient carrying case
- Handy fold-up cart
- Detachable accessory bag



Airline policies for traveling with a portable oxygen concentrator vary. Patients should check with their airlines before flying with SimplyGo.

| Specifications           |  |
|--------------------------|--|
| Oxygen concentrations    | 87 to 96% at all settings  |
| Pulse mode bolus size    | 1=12 ml, 1.5=18 ml, 2=24 ml, 2.5=30 ml, 3=36 ml, 3.5=42 ml, 4=48 ml, 4.5=54 ml, 5=60 ml, 5.5=66 ml, 6=72 ml                            |
| Continuous flow settings | 0.5 LPM, 1 LPM, 1.5 LPM, 2 LPM   |
| Sound level              | ≤ 43 dB at pulse setting 2   |
| Size                     | 11.5 in x 10 in x 6 in<br>(29.2 cm x 25.4 cm x 15.2 cm)  |
| Weight                   | 10 lbs (4.5 kg) with battery   |
| Operating temperature    | 41° to 104°F (5° to 40°C)  |
| Storage temperature      | -4° to 140°F (-20° to 60°C) unit only  |
| Altitude                 | Up to 10,000 ft. (3048 m)  |
| Filtration               | Internal bacteria filter   |
| DC power                 | Input 12 to 19 VDC 8.0 A max   |
| AC power                 | Input 100 to 240 VAC, 50/60 Hz, 2.0 A<br>Output 19 VDC 7.9 A max   |
| Power consumption        | 150W while charging,<br>120W while not charging  |
| Alarm types              | Low oxygen purity, No breath,<br>High breath rate, No flow, Low<br>battery, Depleted battery, Wrong<br>battery, System technical fault |
| Lithium Ion battery      | 1.5 lbs (0.7 kg), charge time is<br>approximately 2 to 3 hours per battery<br>from fully discharged                                    |

| Ordering information  |         |
|---|---------|
| <b>Basic system</b>   |         |
| Includes SimplyGo portable oxygen concentrator, carrying case, mobile cart, one rechargeable Lithium Ion battery, DC and AC power cords, accessory bag, and user manual | 1068987 |
| <b>Accessories</b>  |         |
| SimplyGo carrying case  | 1082663 |
| Lithium Ion battery   | 1082662 |
| Battery charger   | 900-103 |
| DC power supply   | 1083692 |
| AC power supply and cord  | 1082661 |
| Airline power cord  | 1083693 |
| Accessory bag   | 1083696 |
| Mobile cart   | 1074885 |
| User manual   | 1083700 |
| Humidifier pouch  | 1083699 |
| SimplyGo humidifier connecting tube kit   | 1101172 |
| SimplyGo humidifier connecting tube kit and pouch   | 1101602 |
| <b>Warranty options</b>   |         |
| SimplyGo two-year total warranty  | 1097623 |
| SimplyGo three-year total warranty  | 1097577 |
| SimplyGo five-year total warranty   | 1097778 |



**Philips Healthcare is part of  
Royal Philips Electronics**

**How to reach us**

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Asia  
+49 7031 463 2254

Europe, Middle East, Africa  
+49 7031 463 2254

Latin America  
+55 11 2125 0744

North America  
+1 425 487 7000  
800 285 5585 (toll free, US only)

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Customer Service  
+1 724 387 4000  
800 345 6443 (toll free, US only)

Philips Respironics International  
Headquarters  
+33 1 47 28 30 82

Philips Respironics Asia Pacific  
+65 6882 5282

Philips Respironics Australia  
+61 (2) 9947 0440  
1300 766 488 (toll free Australia only)

Philips Respironics China  
+86 400 828 6665  
+86 800 828 6665

Philips Respironics Deutschland  
+49 8152 93 06 0

Philips Respironics France  
+33 2 51 89 36 00

Philips Respironics Italy  
+39 039 203 1

Philips Respironics Sweden  
+46 8 120 45 900

Philips Respironics Switzerland  
+41 6 27 45 17 50

Philips Respironics United Kingdom  
+44 800 1300 845

[www.philips.com/respironics](http://www.philips.com/respironics)

Respironics, EverFlo, and SimplyGo are trademarks of Respironics, Inc. and its affiliates. All rights reserved.

Please visit [www.philips.com/simplygo](http://www.philips.com/simplygo)



© 2012 Koninklijke Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Hoehc WMB 11/12/12 MCI 4105293 PN 1097418